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Voice unit module to activate voice communication, activate and adjust the music channel and announcements, equipped with a flat cable for connecting to the display module, complete with single base for surface mounting, white.

The device, installed in the single room and powered directly by the display module 02081.AB, enables hands-free communication between patient and nurse and between nurses; via the voice unit module it is moreover possible to make room, ward and general announcements and broadcast a music channel with the possibility of adjusting the listening volume.

The device is equipped with 4 front buttons to activate voice communication, switch on, off and adjust the volume (decrease and increase) of the music channel. It is connected to the display module 02081.AB by means of the flat cable supplied.

CHARACTERISTICS.

- Rated supply voltage (from display module 02081): 5 V d.c. \pm 5%.
- Absorption: 5 mA.
- Speaker output power: 0.15 W/16 Ω .
- Speakers: 2 x 8 Ω -250 mW in series.
- Operating temperature: +5 $^{\circ}$ C - +40 $^{\circ}$ C (indoor).

INSTALLATION.


Vertical installation with double base:

- to carry out semi-recessed mounting on light walls, on boxes with 60 mm distance between centres or on 3-gang boxes, use the double base;
- connect the flat cable to the voice unit module 02082.AB and hook it onto the double base 02083 taking care to lay out the cable;
- before hooking the display module 02081.AB onto the double base 02083, connect the extractable terminals 02085 (bus, inputs/outputs, + OUT and -).

Vertical/horizontal installation with single base:

- use the single base to perform the installation;
- connect the flat cable to the voice unit module 02082.AB and hook it onto the single base taking care to lay out the cable;
- before hooking the display module 02081.AB onto its single base, connect the extractable terminals 02085 (bus, inputs/outputs, + OUT and -).

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The voice unit module is used to perform the following functions:

Voice communication.

Systems equipped with speaker modules enable remote communication between the rooms which have been provided with nurse present signalling (green button on the display module) or between the supervisor and the room which has been provided with presence signalling. The volume level of the voice unit cannot be changed from the terminal.

- Pressing button H  once only (fully illuminated) starts hands-free communication with the terminal from which the call is made;
- on pressing button H  a second time (minimum illumination) the voice communication is interrupted.

If there is more than one call, with button A  of the display module 02081.AB, it is possible to scroll through the list of these calls and select the one you want to answer.

Button E  lights up fully when a call is made to the room (for example via VOX) or when there is voice communication; in the case of voice communication piloted by the nurse,  lights up to indicate you can speak (voice unit module in transmission).

The "direction" in which the communication is made is indicated by the selfsame button (button E  on = speak; button E  off = listen).

The mode with which this communication is managed (full duplex/half duplex) is established by the device that triggers it:

- telephone coupler ⇒ always full duplex;
- voice ⇒ depending on the chosen configuration.

In the latter mode, the half duplex switching can take place in two ways:

- **Hands-free**, where the "direction" of the communication is established by the tone of voice; the exchange is made when the voice unit module recognizes a higher sound level of one speaker rather than the other. This type of solution is used in rooms that are not very noisy.

- **Push to talk**, where the communication exchange between the speakers takes place by pressing the button E  (press to speak, release to listen) by the health personnel present in the control room or in a room where assistance is being given; the switching is controlled by the terminal that requested the voice unit connection. This type of application is used in noisy rooms.

Music transmission.

Connecting the system to an audio source, when the system contains a phone coupler, the voice unit modules enable transmitting a music channel.

- Pressing button E  switches the music transmission on and off (the button is illuminated);
- pressing button F  lowers the volume;
- pressing button G  increases the volume.

Buttons E , F , G , and H  are backlit with a red light for location in the dark.

When the voice or music channel is activated, the display will show the symbol  with the set volume level .

INSTALLATION RULES.

Installation should be carried out by qualified staff in compliance with the current regulations regarding the installation of electrical equipment in the country where the products are installed.

Recommended installation height: from 1.5 m to 1.7 m.

CONFORMITY.

EMC directive.

Standards EN 60950-1, EN 61000-6-1, EN 61000-6-3.

REACH (EU) Regulation no. 1907/2006 – Art.33. The product may contain traces of lead.

WEEE - Information for users

 If the crossed-out bin symbol appears on the equipment or packaging, this means the product must not be included with other general waste at the end of its working life. The user must take the worn product to a sorted waste center, or return it to the retailer when purchasing a new one. Products for disposal can be consigned free of charge (without any new purchase obligation) to retailers with a sales area of at least 400 m², if they measure less than 25 cm. An efficient sorted waste collection for the environmentally friendly disposal of the used device, or its subsequent recycling, helps avoid the potential negative effects on the environment and people's health, and encourages the re-use and/or recycling of the construction materials.